EV, Rêma | Patient Transfer Arm Mobile Rocker Recliners - 24" Seat Width - 550 lb Weight Rating







Rêma | R1127L | RC1127L 32.5" Overall Width

DIMENSIONS | E1127L | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	25
Depth	38	Seat Depth	21.5	Arm Height	25
Height	44.75	Seat Height	19.75	Weight	151
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	20.5

DIMENSIONS | R1127L | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	23.75
Depth	38	Seat Depth	21.5	Arm Height	25
Height	43.5	Seat Height	19.75	Weight	153
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	19.5


1127L MOBILE ROCKER RECLINER INFINITE POSITION LOCKING BACK - LEFT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM 4 5 6 7 8 9 10 11 4019 4169 4219 4379 4529 4689 4939 5219 5629


C1127L MOBILE CHAISE ROCKER RECLINER INFINITE POSITION LOCKING BACK - LEFT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11
4269	4479	4549	4689	4869	5089	5299	5629	6029

Option	S	List Price	Option	S	List Price
U W J	Urethane Arm Caps* Wood Arm Caps* Solid Surface Arm Caps* (Not Available on Leda)	Standard \$199 \$499	MBB MBS FBH 3S4	Moisture Barrier Back Moisture Barrier Seat Foley Bag Holder All 3" Swivel Locking Caster	\$50 \$50 \$79 \$19
Z CAL	Multi Cover CAL TB133	Use Highest Grade See page 12	AHP	Adjustable Headrest Pillow (*Requires 1 yd of fabric)	\$299



Infinite Position Locking Back Lever (1127 + 4127 Models Only) Infinite Position Locking Back lever allows patient/caregivers to lock back in any recline position. Locking Back Lever is located on the inside of the left arm (facing) on select models. Locking Back Lever is located on the outside of the left arm (facing) on Dual Lateral Transfer Arm models.



Rocker/Lock Mechanism (1127 + 4127 Models Only)
The mobile rocker recliner features a locking mechanism that locks the rocker making it easy on the patient for entering/exiting the chair and also during transport.

Please refer to pages 160/161 for Options descriptions & images. Please refer to page 158 for the recliner ordering guide.







^{*}Please refer to page 14 for finishes.