EV, Rêma | Patient Transfer Arm Mobile Rocker Recliners - 24" Seat Width - 550 lb Weight Rating







Rêma | R1127R | RC1127R 32.5" Overall Width

DIMENSIONS | E1127R | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	25
Depth	38	Seat Depth	21.5	Arm Height	25
Height	44.75	Seat Height	19.75	Weight	151
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	20.5

DIMENSIONS | R1127R | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	23.75
Depth	38	Seat Depth	21.5	Arm Height	25
Height	43.5	Seat Height	19.75	Weight	153
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	19.5

1127R MOBILE ROCKER RECLINER INFINITE POSITION LOCKING BACK - RIGHT FACING FOLD DOWN ARM - URETHANE ARM CAPS

7 COM 5 6 8 9 10 4 11 4939 4019 4169 4219 4379 4529 4689 5219 5629

C1127R MOBILE CHAISE ROCKER RECLINER INFINITE POSITION LOCKING BACK - RIGHT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11
4269	4479	4549	4689	4869	5089	5299	5629	6029

Options		List Price	Options	S	List Price
U W J	Urethane Arm Caps* Wood Arm Caps* Solid Surface Arm Caps* (Not Available on Leda)	Standard \$199 \$499	MBB Moisture Barrier Back MBS Moisture Barrier Back FBH Foley Bag Holder 3S4 All 3" Swivel Locking Caster AHP Adjustable Headrest Pillow (*Requires 1 yd of fabric)	Moisture Barrier Back	\$50 \$50 \$79 \$19
Z CAL	Multi Cover CAL TB133	Use Highest Grade See page 12		,	\$299



Infinite Position Locking Back Lever (1127 + 4127 Models Only) Infinite Position Locking Back lever allows patient/caregivers to lock back in any recline position. Locking Back Lever is located on the inside of the left arm (facing) on select models. Locking Back Lever is located on the outside of the left arm (facing) on Dual Lateral Transfer Arm models.



Rocker/Lock Mechanism (1127 + 4127 Models Only)
The mobile rocker recliner features a locking mechanism that locks the rocker making it easy on the patient for entering/exiting the chair and also during transport.

*Please refer to page 14 for finishes.

Please refer to pages 160/161 for Options descriptions & images. Please refer to page 158 for the recliner ordering guide.





