





## EV, Rêma, Leda, Koncert | Mobile Rocker Recliners - 24" Seat Width - 550 lb Weight Rating

EV | E1127 | EC1127  
32.5" Overall WidthRêma | R1127 | RC1127  
32.5" Overall WidthLeda | L1127 | LC1127  
32.5" Overall WidthKoncert | K1127 | KC1127  
32.5" Overall Width

## DIMENSIONS | R1127 | L1127 | K1127 | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	23.75
Depth	38	Seat Depth	21.5	Arm Height	25
Height	43.5	Seat Height	19.75	Weight	153
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	19.5

## DIMENSIONS | E1127 | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	25
Depth	38	Seat Depth	21.5	Arm Height	25
Height	44.75	Seat Height	19.75	Weight	151
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	20.5

## 1127 MOBILE ROCKER RECLINER INFINITE POSITION LOCKING BACK - UPH ARMS

COM	4	5	6	7	8	9	10	11
3369	3519	3669	3859	4089	4359	4619	4969	5689

## C1127 MOBILE CHAISE ROCKER RECLINER INFINITE POSITION LOCKING BACK - UPH ARMS

COM	4	5	6	7	8	9	10	11
3599	3759	3899	4099	4329	4589	4859	5209	5919

## Options

F	Upholstered Arms
U	Urethane Arm Caps*
W	Wood Arm Caps*
J	Solid Surface Arm Caps*
	(Not Available on Leda)
Z	Multi Cover
CAL	CAL TB133


## List Price

Standard
\$149
\$199
\$499
Use Highest Grade
See page 12

## Options

MBB	Moisture Barrier Back	\$50
MBS	Moisture Barrier Back	\$50
FBH	Foley Bag Holder	\$79
3S4	All 3" Swivel Locking Caster	\$19
AHP	Adjustable Headrest Pillow	\$299
	(*Requires 1 yd of fabric)	


## List Price



## Infinite Position Locking Back Lever

(1127 + 4127 Models Only)

Infinite Position Locking Back lever allows patient/caregivers to lock back in any recline position. Locking Back Lever is located on the inside of the left arm (facing) on select models. Locking Back Lever is located on the outside of the left arm (facing) on Dual Lateral Transfer Arm models.

Rocker/Lock Mechanism  
(1127 + 4127 Models Only)

The mobile rocker recliner features a locking mechanism that locks the rocker making it easy on the patient for entering/exiting the chair and also during transport.

\*Please refer to page 14 for finishes.

Please refer to pages 160/161 for Options descriptions & images.  
Please refer to page 158 for the recliner ordering guide.

